

	Document Owner Department: QC	BT-SPEC-0060
		Page 1 of 5
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
BRILLIANT GREEN BILE 2% BROTH CM0031		

BRILLIANT GREEN BILE 2% BROTH

CM0031

Formula

Peptone	grams per litre	10.0
Lactose		10.0
Ox-Bile (purified)		20.0
Brilliant green		0.0133

Directions

Add 40g to 1 litre of distilled water. Mix well and distribute into containers fitted with Durham's tubes. Sterilize by autoclaving at 121°C for 15 minutes. Double strength broth - heat the dissolved broth at 100°C for 30 minutes - do not autoclave.

Physical Characteristics

Pale green, free flowing powder
 Colour on reconstitution - green
 Moisture level less than 7%
 pH - 7.4 ± 0.2 at 25°C
 Clarity - clear

Microbiological Tests Using Optimum Inoculum Dilution

Control Media : Tryptone Soya Agar or Columbia Blood Agar Base enriched with 5% v/v horse blood, where appropriate.

ThermoFisher S C I E N T I F I C	Document Owner Department: QC	BT-SPEC-0060 Page 2 of 5
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
BRILLIANT GREEN BILE 2% BROTH CM0031		

Reactions after incubation at 30 ± 2°C for 24 ± 2 hours

Medium is challenged with 10-100 colony forming units.

<i>Enterobacter aerogenes</i>	NCTC 9735	Turbid growth and gas
-------------------------------	-----------	-----------------------

A satisfactory result is represented by visible growth and gas.

Reactions after incubation at 30 ± 2°C for 48 ± 2 hours

Medium is challenged with 1E+04 to 1E+06 colony forming units.

<i>Staphylococcus aureus</i>	ATCC® 25923	No growth
------------------------------	-------------	-----------

Negative strains are inhibited.

Reactions after incubation at 37 ± 2°C for 24 ± 2 hours

Medium is challenged with 10-100 colony forming units.

<i>Escherichia coli</i>	ATCC® 8739	Turbid growth and gas
<i>Escherichia coli</i>	ATCC® 25922	Turbid growth and gas
<i>Enterobacter aerogenes</i>	NCTC 9735	Turbid growth and gas
<i>Citrobacter freundii</i>	ATCC® 43864	Turbid growth and gas

A satisfactory result is represented by visible growth and gas.

Reactions after incubation at 37 ± 2°C for 24 ± 2 hours

Medium is tested at double strength

Medium is challenged with 10-100 colony forming units.

<i>Escherichia coli</i>	ATCC® 8739	Turbid growth and gas
-------------------------	------------	-----------------------

A satisfactory result is represented by visible growth and gas

ThermoFisher SCIENTIFIC	Document Owner Department: QC	BT-SPEC-0060 Page 3 of 5
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
BRILLIANT GREEN BILE 2% BROTH CM0031		

Reactions after incubation at 37 ± 2°C for 48 ± 2 hours

Medium is challenged with 1E+04 to 1E+06 colony forming units.

<i>Staphylococcus aureus</i>	ATCC® 25923	No growth
<i>Enterococcus faecalis</i>	ATCC® 29212	No growth or turbid growth no gas
<i>Enterococcus faecalis</i>	ATCC® 19433	No growth or turbid growth no gas

E. coli confirmation test

Reactions after incubation at 44 ± 1°C for 24 ± 2 hours

<i>Escherichia coli</i>	ATCC® 25922	Turbid growth and gas
<i>Escherichia coli</i>	ATCC® 8739	Turbid growth and gas
<i>Enterobacter aerogenes</i>	NCTC 9735	No growth or turbid growth no gas
<i>Citrobacter freundii</i>	ATCC® 43864	No growth

Gram +ve sporing anaerobes test

Reactions after incubation at 30 ± 2°C for 48 ± 2 hours

Enriched with 10% v/v pasteurised milk

<i>Clostridium perfringens</i>	ATCC® 13124	No gas
--------------------------------	-------------	--------

Reactions after incubation at 37 ± 2°C for 48 ± 2 hours

<i>Clostridium perfringens</i>	ATCC® 13124	No gas
--------------------------------	-------------	--------

Reactions after incubation at 44 ± 2°C for 48 ± 2 hours

<i>Clostridium perfringens</i>	ATCC® 13124	No gas
--------------------------------	-------------	--------

ThermoFisher S C I E N T I F I C	Document Owner Department: QC	BT-SPEC-0060 Page 4 of 5
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
BRILLIANT GREEN BILE 2% BROTH CM0031		

Testing performed in accordance with ISO11133:2014

Reactions after incubation at 30 ± 2°C for 24 ± 2 hours

Medium is challenged with 10-100 colony forming units

<i>Escherichia coli</i>	ATCC® 8739	WDCM00012	Turbid growth and gas
<i>Escherichia coli</i>	ATCC® 25922	WDCM00013	Turbid growth and gas
<i>Citrobacter freundii</i>	ATCC® 43864	WDCM00006	Turbid growth and gas

A satisfactory result is represented by visible growth and gas.

Reactions after incubation at 30 ± 2°C for 48 ± 2 hours

Medium is challenged with 10-100 colony forming units

<i>Enterococcus faecalis</i>	ATCC® 29212	WDCM00087	No growth to turbid growth no gas
<i>Enterococcus faecalis</i>	ATCC® 19433	WDCM00009	No growth to turbid growth no gas

Negative strains are inhibited or shall produce partial inhibition and no gas.

	Document Owner Department: QC	BT-SPEC-0060
		Page 5 of 5
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
BRILLIANT GREEN BILE 2% BROTH CM0031		

Revision History

Section / Step	Description of Change	Reason for Change	Reference
Entire Document	Update to test specification	Change control	BT-CC-1506