



Tubes, culture

PYREX®

- Ideal for culture work
- Manufactured from Pyrex® borosilicate glass
- High resistance to attack from water reduces leaching of contaminants which can cause pH changes
- Phenolic screwcaps with inert PTFE lined rubber disc
- Cap can withstand repeated dry heat and steam sterilisation cycles (121°C for 2 hours)

Product Code	O.D. x length mm	Wall thickness mm	Thread size mm	Pack qty.
1636/24MP	14 x 100	1.5	13	40
1636/26MP	16 x 100	1.8	15	40
1636/04MP	16 x 125	1.8	15	40
1636/30MP	16 x 160	1.8	15	40
1636/32MP	18 x 100	1.8	18	40
1636/34MP	18 x 180	1.8	18	40
1636/36MP	20 x 100	1.8	20	40
1636/38MP	20 x 200	1.8	20	40
1636/10MP	22 x 200	1.8	22	40
1636/42MP	26 x 100	2.0	25	40
1636/44MP	26 x 200	2.0	25	40



*Culture tubes are not recommended for centrifuge work.
Do not autoclave with a tightly screwed cap, always slacken off. Failure to adopt this procedure may result in tube breakage.*



Tubes, culture, spares

Product Code	Description	Thread size mm	Pack qty.
4503/13	Cap with insert	13	20
4503/15	Cap with insert	15	20
4503/18	Cap with insert	18	20
4503/20	Cap with insert	20	20
4503/22	Cap with insert	22	20
4503/25	Cap with insert	25	20
4513/13M	Spare insert	13	20
4513/15M	Spare insert	15	20
4513/18M	Spare insert	18	20
4513/20M	Spare insert	20	20
4513/22M	Spare insert	22	20
4513/25M	Spare insert	25	20