

Cleaning

Tissues, Wipes, Towels and Cloths continued

All dimensions and quantities are nominal.

Towels, Kleenex, folded

Highly absorbent, 3-ply "airflex" fabric which has excellent wet strength and softness. Towel size 220 x 320mm. Supplied in self-presenting packs of 96 towels.

CL730-26 Folded towels, 6771

CL732-10 Dispenser, plastic, wall mounting

Wipes, cleanroom, Kimtech

Precision wipes. Suitable for use in ISO cleanroom Class 7 environments and above. Strong, low-linting, absorbent and chemically inert. Resistant to solvents and most chemicals. Supplied in wipe sizes and packs as indicated.

	Ref.	Wipe, mm	Pack
CL734-15	7623	340 x 380mm	roll 600
CL734-25	7624	350 x 385mm	35

Duster, yellow

Soft with close texture. Dimensions 500 x 350mm. Supplied in pack of 10.

CL742-10 Dusters

Polishing Cloths, Selvyt®

The original soft texture polishing cloth which can be washed repeatedly. Supplied in packs of 10.

CL744-11 250 x 250mm

CL744-14 350 x 350mm

Wiping Cloths, J-type

Bonded cellulose fibre. Can be used for wet or dry cleaning. Dimensions 330 x 500mm. Supplied in pack of 50.

CL745-40 Wiping cloths

Glass Cloth, Linen

For drying of glassware and general equipment. Dimensions 750 x 500mm. Supplied in pack of 10.

CL745-60 Glass cloths

Wipes, bactericidal, Sani-Cloth 70

For disinfection of surfaces in clinical environments. Perforated roll of non-woven material saturated with 70% v/v Isopropyl alcohol B.P. Wipe dimensions 200 x 220mm. Supplied in a dispenser tub of 200. **Please note:** these items are classed as hazardous material for export and require packing as such.

CL745-65 Wipes, bactericidal

Wipes, bactericidal, Sani-Cloth 70, sterile

Generally as CL745-65 but supplied individually wrapped in sachets and sterile. Wipe dimensions 200 x 128mm. Supplied in a pack of 100 wipes.

CL745-70 Wipes, bactericidal, individually wrapped, sterile

Lens Tissues

Tissues, Whatman 105

For lens cleaning.

	Dimensions, mm	Pack qty.
CL760-10	100 x 150	625
CL760-15	200 x 300	100
CL760-20	460 x 570	500



CL745-70



CL730-26 with CL732-10



CL734-25



CL745-65



CL734-15



CL760